

P: +1.206.365.1256
 E: cservice@usbiotek.com
 1620 Linden Av N
 Shoreline WA, 98133
 Dr.- US BIOTEK

LAB ID : 3888853
 Collection Date : 08-May-2023
 Received Date: 08-May-2023

Accession #: 0000000000

GIT ASSESSMENTS

STOOL, SPOT

Result	Range	Units
>500	> 200	ug/g



Pancreatic Elastase 1

Digestive Markers Comment

PANCREATIC ELASTASE: Normal exocrine pancreatic function.
 Pancreatic Elastase reflects trypsin, chymotrypsin, amylase and lipase activity.
 This test is not affected by supplements of pancreatic enzymes.
 Healthy individuals should be producing >500 ug/g of PE-1 under normal/healthy conditions.
 PE-1 levels between 200 - 500 ug/g may indicate suboptimal production.
 PE-1 levels <200 ug/g indicate clear inadequate production.
 The clinician should therefore consider digestive enzyme supplementation if one or more of the following conditions is present:
 Loose watery stools, Undigested food in the stools, Post-prandial abdominal pain, Nausea or colicky abdominal pain, Gastroesophageal reflux symptoms, Bloating or food intolerance.

Methodology: EIA

Tests ordered: 2014

DISCLAIMER:
 Any information provided by us is for information purposes only.
 Commentary is provided to the practitioner for educational purposes and should not be interpreted as diagnostic or as treatment recommendations.
 Diagnosis and treatment decisions are the practitioner's responsibility.